

3M Science.
Applied to Life.™

3M™ Scotchbond™ Universal Plus
Adhesive

**More than just
a strong bond.**



3M™
Scotchbond™
Universal Plus
Adhesive takes
universal adhesives to
the next level.
It keeps all the benefits of the
original 3M™ Scotchbond™
Universal Adhesive, the first truly
universal adhesive, while adding
even greater control and
predictability – because a good
adhesive should provide more
than just a strong bond.

From excellent...

Built on the clinically proven and world's most researched universal adhesive.

-  **All direct and indirect indications**
light-curing, self-curing or dual-curing
methacrylate-based composites / cements
-  **All etching techniques**
total-etch, self-etch, and selective
enamel etch
-  **Bonds to all restorative materials**
without the need for separate primers
-  **Available in vial and unit dose**
for efficient hygiene management
-  **Virtually no post-operative sensitivity**
in field evaluation:

99% of dentists reported
that **3M™ Scotchbond™ Universal Plus**
Adhesive shows virtually no post-
operative sensitivity.*

*Field evaluation with 309 dentists and over 20,000 applications

...to exceptional


Formulation carefully optimised to achieve unique features.

-  **First radiopaque universal adhesive**
reducing the risk of X-ray misdiagnosis and
invasive overtreatment*
-  **Bonds and seals caries-affected dentine**
maximising preservation of natural tooth structure,
supporting minimally invasive preparation
guidelines*
-  **Gold standard adhesion to all dental substrates**
including glass ceramics
-  **Full dual- and self-cure compatibility**
no need for an additional dual-cure activator bottle
-  **BPA free formulation**
free of BPA derivatives like
BisGMA
-  **Fully aligned system with 3M™ RelyX™
Universal Resin Cement**
excellent bond strength for virtually all resin
cement indications.* Adhesive cured by cement.



Added radiopacity


On average, dentists see at least one hard to interpret X-ray per day due to radiolucency under an existing restoration. This could lead to misdiagnosis and invasive overtreatment. 3M™ Scotchbond™ Universal Plus Adhesive is the first radiopaque universal adhesive with dentine-like radiopacity, helping to overcome this issue.*



Non-Radiopaque Adhesive

3M™ Scotchbond™ Universal Plus Adhesive

Intentional adhesive pooling under Class I restorations – X-rays of extracted human molars, filling material 3M™ Filtek™ One Bulk Fill Restorative. 3M internal data, H. Loll et al., J Dent Res 99 (Spec Iss A): No. 757, 2020



“Radiolucent adhesive pooling underneath a restoration doubles the risk of restoration replacement.”

**Dr. John O. Burgess,
University of Alabama, USA**

Ordering information

Item no. Product information

3M™ Scotchbond™ Universal Plus Adhesive – Vial

| | |
|--------------|---|
| 41293 | Intro Kit Vial 1 Vial 5ml, 50 Disposable Brushes, 1 Etchant Syringe 3ml, 25 Dispensing Tips |
| 41294 | Refill Vial – 1 × 5ml |
| 41295 | Refill Vial – 3 × 5ml |
| 41263 | 3M™ Scotchbond™ Universal Etchant Etching Gel 2 Syringes x 3ml, 50 Dispensing Tips |



Vial
with smooth and even surface –
easy to wipe off excess

3M™ Scotchbond™ Universal Plus Adhesive – Unit Dose

| | |
|--------------|--|
| 41297 | Intro Kit Unit Dose 50 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips |
| 41298 | Unit Dose 100 100 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips |
| 41299 | Unit Dose 200 200 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips |
| 41304 | Unit Dose 400 400 unit doses, 1 Etchant Syringe 3ml, 25 Dispensing Tips |



Unit Dose
for efficient hygiene management –
reduces the risk of cross
contamination

*3M internal data



3M Oral Care
3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough
Leicestershire, LE11 5RB
+0800 626 578

3M Ireland Ltd
The Iveagh Building
Carrickmines Park
Carrickmines Dublin 18
Ireland
+353 1 280 3555



3M.co.uk/Scotchbond-plus

☎ 08705 360036

✉ OralCareUK@mmm.com

📱 @3MOralCareUK

3M, Filtek, RelyX and Scotchbond are trademarks of 3M Company. © 3M 2020. All rights reserved. G0113A / 17422317_UK20201111